



D3050 Root Tip Extraction Forceps: Used for removing root tips that may have been broken off during tooth extraction & are inside the tooth socket where regular extraction forceps cannot access. They are also useful for small teeth when extracting such as in felines. They are not general purpose extraction forceps for larger teeth when the regular & larger extraction forceps are required & should be used.



D3060 Fragment Extraction Forceps: Ideal for picking up small fragments of tooth, root, bone etc that may be loose & not still attached & floating in a small tooth socket & cannot be accessed by either Root tip & or regular sized extraction forceps. They may be used for the extraction of the smallest of delicate teeth.



D10100 iM3 Explore/Measuring Probe Dog:




Explorer/Measuring Probe
The probe end has periodontal markings etched in black. The fine but robust tip is ideal for periodontal probing in cats or dogs. The markings will not fade even after many autoclave cycles. The other end of the probe has a shepherd hook explorer to detect tooth defects both supra and subgingivally.

Cat. No
D1010

D1010 Explorer/Measuring Probe Order Code: D10100

The probe end has periodontal markings 1,2,3,5,7,8,9,10 mm printed on the handle. The fine but robust tip is ideal for periodontal probing in the cat or dog. The markings are indented into the probe and will not fade even after many autoclave cycles. The other end has a shepherd hook explorer to detect tooth defects both supra and subgingivally.



● iM3's Latest Offering

iM3's new stainless steel scalpel handle was expressly developed for creating gingival flaps using a #15 blade during dental surgical extractions. The round pen like handle allows for precise control whilst making the incisions and the handle sits comfortably in your closed hand. The new ergonomic design is a great improvement on the traditional flat scalpel handles that are difficult to control.

Order Code: T1400N

