

Dr. Larson's

Aural Haematoma Tube

Primary treatment of aural haematomas is a routine surgical procedure, however, despite the various methods of surgical treatments developed, recurrence is common. Treatment methods involving continued drainage of haematoma fluid production have been developed. The Dr. Larson's™ Aural Haematoma Tube is a preferred method for aural haematoma treatment.

To drain the haematoma, surgically insert the Dr. Larson's™ Aural Haematoma Tube at the most distal or gravity dependent part of the haematoma (lowest point when the ear flap is hanging down). Leave the tube in place for 10-14 days. The retaining fingers will secure the tube in place during the healing and treatment period.



- Ideal design for use as a drain tube for auricular haematomas in dogs and cats
- Flexible retaining fingers secure the tube in place without the use of tape
- Flow-through design of tube allows for continuous drainage of affected area and infusion of medicine or flushing solution
- Screw-on cap may be easily removed for draining or medicating the haematoma
- Individually packaged and gamma irradiated

Tube placed at distal part of haematoma



Item	Description	UOM
AHT	Dr. Larson's™ Aural Haematoma Tube	Bag of 10

Distributed by:

ZebraVet

Ph: 1800 477 385

www.zebravet.com.au