

Panania Vet Clinic, NSW - Leanne, Practice Manager
Not pulling apart at connection for portals. 2 portals as compared to some other giving set.



génia AUSTRALIA

KEEP CALM AND LOVE YOUR VET

IV SETS

<p>Injection Ports</p>	<p>Luer-lock</p> <p>Luer-lock connection</p>	<p>Clamp</p> <p>Flow Regulator</p>		
<p>PERFUVET® DUO Premium giving set Calibration: 20 drops/ml</p>	<p>PERFUVET® PLUG 2 in 1: admin set + 25 cm extension Calibration: 20 drops/ml</p>	<p>PERFUVET® MICRO IV set with flow regulator Calibration: 60 drops/ml <i>Ideal for paediatrics or when running out of pump.</i></p>	<p>INFUJECT-SITE Extensions sets: 26.5 cm or 100 cm</p>	<p>PERFUVET® TWIST Tangle free IV set device. Spiral tubing avoids kinking.</p>

2 Injection Ports: placed at 10 cm and 110 cm from the catheters.
 Compatibles with most IV pumps.

Contact us for more information: contact@geniaaustralia.com - www.geniaaustralia.com