## PERFUSION -

## ■ PERFUVET® TWIST

NEW

Anchored, rotating infusion fixation set that prevents the animal from tangling the tubing or blocking fluid flow. Ideal for fractious hospitalised patients.

- One base to attach to the cage equipped with double-sided tape for secure, fast attachment. Contains four holes to allow screwing into place.
- Two clasps that attach the infusion tubing to the base, allowing unlimited rotation.
- One 3 m length of sterile spiral PVC tubing with a Luer-Lock tip and a ring for attaching to the collar as needed.
- Perfuvet<sup>®</sup> Twist is easily attached to the cage; the base can be left on the cage and reused as often as necessary.
- The unlimited rotation prevents the animal from becoming tangled in the infusion set and blocking flow or requiring a time-consuming change of tubing.
- The spiral tubing allows ease of movement and prevents excessive pulling on the catheter.
- Saves time, improves safety of hospitalised animals, easy to use.
   Use on any size of animal (dogs, cats, rabbits...).

Sold as a kit (One fixed base, one clasp and three sets of single-use spiral tubing).


 Code
 Model
 Pack.

 157450
 Perfuvet<sup>®</sup> Twist kit (I base, 2 clasps and 3 sets of single-use spiral fubring)
 Unit

 157451
 Single-use spiral fubring
 x 3

GONE ARE THE DAYS OF ANIMALS
GETTING TANGLED IN THEIR INFUSION
SET, PERFUVET® TWIST, THE REVOLUTION
IN PET HOSPITALISATION!





Watch the video version on our Youtube channel!



A TIME SAFETY IN SAVER EASY TO USE HOSPITALISATION







INJECTION

