




BleedSolv™ Haemostatic Gauze

BleedSolv™, a TGA approved haemostatic gauze to control bleeding from open wounds and in body cavities. BleedSolv™ is a collagen-like natural substance created from cellulose. On contact with blood, BleedSolv™ expands to 3-4 times its original size, absorbs exudates (fluids) at the wound site and converts to a gel. Because of BleedSolv™'s purity it degrades simply into polysaccharides and water. No need to debride.

A haemostatic gauze applied directly to the site of bleeding to stop bleeding fast in cuts, scratches, burns, incisions, trauma and wounds of all kinds



First wherever bleeding occurs

BleedSolv™ forms a critical part of effective first aid and procedure management. Available in a range of products and sizes for any situation. BleedSolv™ conforms easily to any wound shape and location, is safe for use on the head and in the chest cavity, and is safe for use on animals. BleedSolv™ is particularly beneficial in remote locations distant from advanced care, high-risk situations and getting animals moving earlier.



Effective Care

BleedSolv™ is effective in front-line trauma and first aid care, as well as for dentistry, obstetrical and gynaecological procedures, veterinary procedures and animal husbandry. BleedSolv™ is compact, light and does not require special storage or training. Ideally suited to carry on your person or in mobile situations.




Technology saving lives

BleedSolv™ gauze conforms to the shape of any wound, effective in deep wounds where there is profuse bleeding and disseminated haemorrhage. Following clotting, bioabsorbable BleedSolv™ is rinsed from the wound with saline or water. There is no disturbance of the clotting site to debride the product from the wound. BleedSolv™ accelerates recovery and healing. Essential for remote locations where access to secondary care may be delayed.



FA Ready Pty Ltd



Ahead of time

Sales & Service

Phone: (08) 6143-4600

Mobile: 0415 359 968

Email: service@bleedsolv.com

BleedSolv™

Stop bleeding fast

Topical Haemostatic Dressings

www.bleedsolv.com

BleedSolv™ Haemostatic Gauze


| Product Code | Dimensions | # Items Per Box |
|--------------|--------------|-----------------|
| FAHG-V55 | 5 x 5cm | 12 |
| FAHG-V51 | 5 x 10cm | 12 |
| FAHG-V11 | 10 x 10cm | 12 |
| FAHG-V2020 | 20 x 20cm | 30 |
| FAHG-V75900 | 75 x 900cm | 5 |
| FANG-V2A | Tampon | 10 |
| FADG-V7575 | Dental Gauze | 20 |

NEW **Pet kits now available for retail (FAPK00, FAPK10, FAPK20)



Stop Bleeding Fast

Ahead of time



FA Ready Pty Ltd