


iM3 Airgid – Testimonials

“These newly designed haemostatic sponges are the perfect size for equine dental use and over the past three months we have used them in a number of cheek tooth, incisor and canine alveoli. Placing these immediately post extraction now forms part of our standard procedure for alveolar socket management and so far, we have received excellent feedback from referral clinicians re-examining sockets at a later date”.

Chris Pearce, Equine Dental Specialist, and director at EDC Ltd



<https://www.paperturn-view.com/im3-dental/im3-catalogue-australia?pid=MTA102136&p=47&v=17.1>

iM3 Airgid – Testimonials



“I have used the 20mm³ Airgid sponges on several occasions to pack the alveolus following tooth extraction. I am very pleased with the degree of granulation and lack of unpleasant odour at follow up, compared to other packing methods. I also like the fact that they are dissolvable, so there is no concern about them being placed right at the apex where packing might be difficult to retrieve. Additionally, their size is very convenient for equine use”.

Sam Baker BVSc CertAVP(ED) MRCVS at Apex Equine Dental Surgery (UK)

<https://www.paperturn-view.com/im3-dental/im3-catalogue-australia?pid=MTA102136&p=47&v=17.1>